10

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- *Custom-made implants*

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Custom made nail for knee arthrodesis

Clou pour arthrodèse de genou lg 700 mm <i>Nail for knee arthrodesis plate</i>	
Ø mm	Référence <i>Reference</i>
10	36.045.00-10
11	36.045.00-11
12	36.045.00-12
13	36.045.00-13

Impacteur <i>Impactor</i>	
Référence <i>Reference</i>	
37.595.00-B-M6	
37.595.00-B-M7	
37.595.00-B-M8	
37.595.00-B-M8	

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The Safe Dual core Universal Nail is an orthopedic cannulated nail, perforated in its length, having at the extremities locking holes in perpendicular planes. The preparation of the nail in the operating includes the wrapping with a Esmach band in order to be filled with polymethylmethacrylate cement (PMMA) impregnated with appropriate antibiotics.

Long. en mm <i>Length</i>	Référence
300	36.046.00-30
320	36.046.00-32
340	36.046.00-34
360	36.046.00-36
380	36.046.00-38
400	36.046.00-40

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If necessary, can be reinforced with a metal rod inserted inside the cement, forming a dual-core cement / metal core that increases resistance to fatigue in more than 200 %. This unique model and multiple sizes implant for all bone segments, is intended to be used as a first indication, to open fractures, bone lengthening procedures over nail, conversion of external fixation to nailing and treatment of osteo-articular infection. It can be used in the femur, tibia, humerus and to do arthrodesis of the knee and ankle.

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


Nuno Craveiro Lopes, M.D.

Director of Orthopedic Department


Garcia de Orta Hospital - Almada, Portugal

site : <https://sites.google.com/site/sothgo1>

mail : nlopes@hgo.min-saude.pt

Cas clinique n°1

Clinical Case N°1

Fracture ouverte de grade II, encloué immédiatement. Le clou s'adapte parfaitement à la courbure du tibia.
This patient had a grade II open fracture and was nailed immediately. You can observe the perfect fit of the nail to the shape of the tibia in the profile.

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The patient at six weeks of follow up, do total load, have no signs of infection and returned to work as a karate teacher. At right, the one year follow-up xr.

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

This patient had a sequelae of an infected total knee replacement. After 10 operations, we proceeded to a knee arthrodesis and tibial lengthening with a Ilizarov apparatus.

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By intolerance, the device was replaced by a SAFE DualCore Universal extra-long nail, reinforced only in the area of the knee and tibial lengthening, leaving some flexibility in the proximal zone. At right, the result after 9 months of follow up.

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Preoperatively:

The proper length of nail is determined preoperatively by radiographs of the opposite bone if it is a case of fracture. For reconstruction procedures the measure must be taken in a calibrated X-ray so that precise measurements can be done. The nail should extend well below the fracture or lesion site. Ideally the nail fills the majority of the medullary canal, and the nail should extend nearly the full bone length. Due to caudal curvature of the femur, caution must be used distally so that the nail is not forced through the cranial cortex. It is important that the screws/bolts are at least 2 cm away from fracture site.

Preparing the nail and medullary canal:

After selecting the proper SAFE nail and prepared it with the DualCore, a standard approach to the bone is used. To open the medullary canal for later nail insertion, the initial guide wire is driven normograde or retrograde. Prepare the medullary canal with reamers to 1,5 to 2 mm larger than the SAFE nail diameter. After reduction is manually attained with correct rotational and axial alignment, the nail is introduced normograde or retrograde until the lesion is secure.

Inserting the nail:

1. *After selecting the proper SAFE Nail, attach the proximal guide (ref 36.049.30), using the Hex Driver tighten the internal screw inside the SAFE Nail.*
2. *Attach the Insertion Tool (ref 36.049.10).*
3. *Insert the SAFE Nail into medullary canal either manually or with a hammer. Place the proximal end of the SAFE Nail below the surface of the joint*
4. *Insert Guide Sleeve (ref 36.050.20) into the proximal hole of the guide.*
5. *Use the trocar to mark bone for a starting point for drill bit.*
6. *Insert Drill Guide into Guide Sleeve (ref 36.050.10).*
7. *Place drill Ø 3 mm bit into Drill Guide.*
8. *Drill hole through near cortex, SAFE Nail hole, and then far cortex.*
9. *STOP and CONFIRM you have drill bit through hole with image intensifier*
10. *Remove drill and drill guide from Guide Sleeve*
11. *Using depth gauge measure the depth of the hole. Make note of the length of screw needed.*
12. *Insert proper screw using screwdriver. The screw and screwdriver can be placed down the guide sleeve.*
13. *Repeat Steps 5-12 for the distal hole of the guide jig*
14. *Introduce the two distal locking screws using the image intensifier and the "free hands" technique.*
15. *Closure of the wound is routine and postoperative radio graphs are taken.*

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- Craveiro Lopes, N., Escalda, C.: SAFE Nail. Interlocking Nail with a Long Acting Antibiotic Releasing Core. An in Vitro Experimental Study and Preliminary Clinical Results. <http://pt.scribd.com/doc/59226554>
- Craveiro Lopes, N., Escalda, C.: SAFE Nail. Use of a Interlocking nail with a long acting antibiotic releasing core in patients with high infection risk. Prospective study in 27 cases. <http://pt.scribd.com/doc/92151128>

Designed and developed by: Nuno Craveiro Lopes, M.D., Director of Orthopedic Department, Garcia de Orta Hospital, Almada, Portugal. <https://sites.google.com/site/sothgo1> - nlopes@hgo.min-saude.pt

Follow-up:

Restricted exercise (leash only, no stairs, etc) is advised until 2 month recheck radiographs are obtained.

After 4-6 months when the fracture has healed, if required the nail can be extracted by removing screws, and attaching the Insertion Tool to dislodge nail from osseous or fibrous ingrowth in the nail holes.

Adverse Effects:

These include breakage of the implant due to excessive activity, prolonged loading upon the device, incomplete healing, excessive force exerted on the implant during insertion, implant migration and/or loosening, metal sensitivity or histological or allergic reaction resulting from implantation of a foreign material, pain or discomfort may occur due to the presence of an implant.

Nerve damage, necrosis of bone or bone resorption, necrosis of tissue or inadequate healing may occur if a incorrect surgical procedure is used.

Warnings:

For safe and effective use of this implant, the surgeon must be thoroughly familiar with the implant and the recommended surgical technique for the device. Device breakage or damage can occur when the implant is subjected to increased loading associated with delayed union, nonunion, or incomplete healing. Improper insertion of the device during implantation can increase the possibility of loosening and migration. The patient must be cautioned, in writing, about the use, limitations, and possible adverse effects of this implant, including the possibility of the device failing as a result of loose fixation and/or loosening and excessive activity before the lesions has healed. The patient must be warned that failure to follow postoperative care instructions can cause the implant and/or treatment to fail.

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
Custom-made shoulder prosthesis and humeral reconstruction

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
Created from frontal and profile telemetric X-rays

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Operative technic : watch the video on our website

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Hinged elbow prosthesis for bone reconstruction

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The progress of chemotherapy justifies surgical techniques which improve the survival quality of those stricken with osteosarcomas and Ewing's sarcomas. If the reconstruction techniques are reliable for adults, for children they run into difficulties because of the problem of future disparity of the limbs's length, particularly the reach of the femur's inferior extremity.

This localisation represents more than forty percent of osteosarcomas and close to twenty percent of Ewing's sarcomas. It is because of these cases that we suggest the growing prosthesis, as soon as the foreseeable disparity, due to the suppression of the growth cartilage of the inferior femur, is equal to or greater than 3 cm.

This prosthesis may be used right away after the tumoral resection, or placed secondarily when the disparity starts to interfere in functionality.

Description

The prosthesis is composed of two parts connected by a hinge. Only a hinged joint can avoid the luxation that lengthening would produce.

The tibial piece is not cemented, so as not to sterilize the growth cartilage of the tibia. The smooth stem is small enough to reduce growth cartilage trauma. When setting, perforation of the growth cartilage should be perfectly centered to limit deviations in the Valgus Varus. The tibial plate is equipped with two small anti-rotation flanges, meant to be included in the trench cut out of the articular cartilage and sub-chondral bone above the conjugation cartilage. The tibial plate rests upon the tibial articular cartilage.

The femoral part created to-size is equipped with a telescopic system that you can operate from the interior of the knee under arthroscopy, thus limiting as much as possible the cutaneous opening and the risk of infection (principle complication of other growing prostheses).

The particularly effective anti-twisting system avoids the releasing of metallic particles when pressure is exerted on the prosthesis.

The tail of the femoral prosthesis, according to need, may be cemented or screwed. It should be solidly set.

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Used for positioning the prosthesis' axle. The axle has a screw-foot so that one may hold it with a screwdriver

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Trial axle

Before cementing the femoral piece, it is prudent to verify the length of the limb at the end of the intervention, and for this, use the trial axle. This prosthesis allows one to make up for an over-cutting error peroperatively, by lengthening the prosthesis. Once you have verified the correct length of the limb, you can cement the femoral implant, and assemble the tibial and femoral fragments with the rings and axle. To avoid the axle coming out of place, it is imperative that you place an anti-extrusion clip with clip applying forceps and verify that the clip is well in place, by turning it within its housing. Closing may then follow; do not hesitate to review the transfer of the medial twin when the musculo-cutaneous covering of the implant seems uncertain.

Lengthening technique

This is performed, if possible, after the end of chemotherapy, when the ailing person has recovered normal immunity defense. The arthroscopy is done according to the usual anaesthesia technique on a knee at 40° flexion. The arthroscopy indicates the inter-condylar indentation of the prosthesis in which is located the head of the growing screw. Inserted through a 3 mm incision, the screw-driver must be pushed all the way into the screw-head's recess before being turned.

In order to begin the lengthening, turn the screw-driver counter-clockwise, if it is a question of a left or right prosthesis. When the lengthening becomes too long, the stretched quadriceps muscle tends to lengthen the knee and diminish the flexion. One must not try to lengthen on a knee bent at less than 30°, because the results will be painful and distinguished by a persistent flexion. Respecting this rule also avoids vasculo-nervous complications observed with other models.


When growth has stopped, it is advised to change the growing prosthesis, whose size is no longer adapted to that of the adolescent, to an adult reconstruction model prosthesis.

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
This gauge permits you to accurately choose the point of entry of the prosthesis' tail and to center it in relation to the conjugation cartilage. Doing the tracing on the pre-operative photograph may help the success of this process.

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It permits execution of the final prosthesis' anti-rotational fin casings, paying attention to the position of the rotation. The ideal being to line up the child's anterior tibial tuberosity with the anterior crown of the tibial plate. Caution : the dimensions of the trial tibia are not the same as the final one. The diameter of the stem is 6 mm.

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