
4

Plaques • *Plates*

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The originality of these implants consists in setting any of these plates without changing the epiphysal reduction of the fracture after having temporarily set it with the pins.

It is very often useful to use this unstable means to temporarily set the separated bone fragments. But once the setting is obtained, the pins can hinder positioning the plate on the surface of the bone. There is therefore a risk of losing the anatomical reconstruction, which has been laboriously obtained.

With this product, the aleatory operative time is omitted ; there is indeed an epiphysal plate, which is similar to the epiphysal part of the permanent plate. This plate is actually a guide (drill guide) with openings that allows a temporary setting of the fracture with the pins. These pins will also be useful in guiding the canulated epiphysal screws.

The holes of the drill guide are actually the same distance apart as the epiphysal screw-holes of the osteosynthesis plate. This drill guide is also equipped with a directional pin, which is located at one of its ends. A pin opening is also there to help line up the instrument on the diaphysal axis.

- *The second stage of the operative scenario consists in removing the drill guide and replacing it with the osteosynthesis plate without removing the pins, since the distance between them is same as the distance between the osteosynthesis plate screwholes.*

Each osteosynthesis plate is also equipped with a special opening located at one of its corners specifically for lining up the pin on the diaphysis.

- *The third stage consists of fixing the diaphysal part to the osteosynthesis plate using the usual cortical screws, then in removing the alignment pin guide on the diaphysis when the plate is well in position.*

On the contrary, the pins that go through the holes of epiphysal part and that maintain the setting are left in position for the moment.

- *The fourth stage consists of using the epiphysal pins as a guide for the canulated drill bit, then as a guide for the epiphysal screws, which are canulated as well.*

The pin guides are removed only when the special epiphysal screws are permanently in position.

Thus, this scenario has at no moment endangered the setting of the focus of the metaphyso-epiphysal fracture, regardless of the level of complexity.

The epiphysal plates have been designed and manufactured for any type of epiphyses concerning the principal long bones. The shape of each plate is adapted to that of the bone surface, onto which it will be applied. Therefore, each model has a right and a left side.

A drill guide, whose shape is similar to the epiphysal part of the permanent implant, is joined to each implant. The canulated epiphysal screws are 6.5 mm in diameter and the pin guides, 1.5 mm. The pins can be oriented in the openings of the drill guide either at 90°, or obliquely up to an angle of 15°, thanks to the unique, removable orientation devices.

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● *Grooved, pre-molded plate, respecting the stape of the epiphysis. The plate fits perfectly onto the projecting edge of the humerus, which simplifies surgery.*

● *In complex fractures, a rigid setting is obtained by placing two screws in the epiphysis and a long oblique screw in the medial column.*

● *No internal protrusion of material threatens the ulnar nerve.*

● *3.5 mm diameter cortical screws.*

Note : *A light adjustment with a jigsaw may be necessary for good adaptation to the humerus. Use traditional instruments in order to take into account different morphologies.*

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Fig. 3

Fig. 4

Osteosynthesis of distal humerus fractures with a pre-molded external plate.

This plate can be placed with the patient in a unilateral reclined position and a median posterior access

With regard to the extensor apparatus, the posture depends on the bone lesions

- *If it is a high supra-capitulum fracture, the extensor apparatus must be respected. The lesion is approached passing on both sides of the triceps.*
- *Extra-articular cutting of the olecranon, if it is a supra-capitulum fracture, or a simple supra- and inner-condylar fracture.*
- *In cases when very good articular control is required, the olecranon will be cut at its base.*

Sight of the ulna nerve is required, but its dissection must be atraumatic, and not mobilized. It will not be transposed at the end of the operation.

After cleaning the joint, perform an evaluation of the articular area lesions and the columnar lesions. The reduction is progressive, under direct inspection, being as anatomical as possible, but without persisting in setting miniscule fragments.

We begin by rebuilding the articular mass, then we join it to the epiphysis. The use of small, self-retaining, bone-holding forceps and fine pins are very useful. These often complicated operatory manoeuvres are highly simplified because we can directly reduce the epiphysis fragments onto the external pre-molded plate.


The setting must be stable in order to allow an early rehabilitation, immediatly after the scarring over of the soft tissue.

Extract of the Dr. P.L. HARDY Conference (Hôpital Ambroise Paré - Boulogne France) presented at the 5th PERFECTO Seminar

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● Plaque en "Y" autosecable *Y-shaped plate*

Long. <i>Length</i> en mm	Nb de trous <i>Holes</i>	Référence
92	2	36.310.16
122	4	36.310.18
152	6	36.310.20
182	8	36.310.22
212	10	36.310.24

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Problems such as union or delayed non-union, new fractures caused by the weakness of the bone structure, fractures across the screw holes due to a concentration of pressure, are overcome by the system of claw-plates.

This system minimises the trauma caused by the retention of the reduced position of the bone fracture and, at the same time, ensures an optimal environment in the area of the fracture so that regeneration of the bone takes place. The creation of this plate, which does not require screws, enables the avoidance of any supplementary damage to the fractured bone and insures and preserves blood circulation all the way to the area of the fracture, a feature which is of utmost importance to the primary regeneration of the bone. The design of the piece is such that during the stabilization of the fracture, the bone experiences a physiological micro-movement. This will aid in regaining normal bone structure (lines of pressure).

This system of claw-plates has been used since 1979 and clinical results show that the healing time of the fracture is reduced 25%. The complication rate is very low, even over a long period.

This direct operative technique for inserting the plate saves time.

Indications

Osteosynthesis of fractures using prosthetic stems.

Advantages

- *Simple operative technique which saves time.*
- *Quick strengthening without complications, without having to observe the effects of the protection of tension osteoporosis under the plaque with X-ray tests.*
- *Optimal circumstances maintained throughout process of regeneration.*

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
- *Physiological movement optimising the restoration of the bone structure, which hinders tension protection.*
- *The application of the plate preserves the blood supply and venous drainage.*
- *The endosteal blood supply is kept intact by avoiding screws and nails.*
- *No concentration of tension around the screw holes.*
- *The time of healing is reduced 25%.*
- *Economic treatment of the patient.*
- *Strengthening : from five weeks.*

Conception

The conception of the clamp plate is such that it is paraskeletal when it is positioned on the bone. Contrary to other circumferential devices, it does not sit flush on the bone or against it, and thus does not oppose the periosteal venous drainage system. The plate covers 2/3 of the circumference of the bone and the groove of its spine adds to its rigidity. The teeth, which have been sharpened and turned toward the interior, come from this part. These teeth are inserted in or around the bone with folding forceps and the the maximum penetration is only one or two millimeters, thus preserving the endosteal vessels.

During the bending stage, the plate takes on an oval form. This creates a space between the plate and the bone, the only point of contact being the "teeth". This minimum contact leaves the periosteum, the muscles and the tendons free.


Responses to Your Questions

Basic philosophic concept

The clamp plate is conceived to inflict minimal trauma when a fracture must be reduced surgically. This plate aims to maintain an optimal healing environment for the fracture, while keeping the reduced position of the operation. Since you do not use a screws to fasten the plaque to the bone, the arterial endosteal blood flow to the creative osteoblastic cells is protected. Moreover, since the plate does not sit flat on the bone, it protects the highly important drainage system via the periosteum up to the tissues.

Fuction of the periosteum

The two most important functions of the periosteum are, first, to aid the bone in draining the blood of the cortex via the periosteum to the tissue, and, second, to aid the healing of the bone at the area of the fracture.

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Stability

The claw-plate stresses the biological aspect of the healing of the bone by preserving some of the more important known factors in the process of healing. Thus, the first objective of the plate is not the mechanic rigidity as one would expect in an apparatus of rigid internal setting with screws. Nevertheless, if the recommended technique is followed, that is to say three pairs of teeth on each side of the fracture line, as well as the correct plaque size for a specific bone, sufficient stability is achieved.

Placing the plate

As the plate changes and passes from a perfect semicircle to a soft oval, it lifts itself slightly in relation to the bone. This small space which rests after the flange plate has been removed is enough to allow the viability of the periosteum. This is proven by an eventual development of new bone, which forms below the plate. This is extra bone, and adds to the rigidity of the whole. This contrasts strongly with a screwplate, in which the bone under the plate becomes vascular, resulting in a deposit due to the blockage of the venous drainage vessels.

Strength of the bone

Since this plaque is not a rigid internal setting device, it may not be able to endure all the constraints on the bone as one sees with a dynamic compression plate with screws. A certain margin of physiological movement or micromovement is granted at the location of the fracture. Even in follow-up cases over ten years, no bone weakening or osteopenial disuse or protection of tension was observed under the flange plate. In fact, the cortex increases in thickness, and because of the supplementary layer of bone that is formed under the spine of the plaque, the whole becomes more solid. Once more, this is the opposite of what one may see with a rigid internal setting device, where the protection of the tension and osteopenia can be noted even at an early stage after the plaque has been set. Moreover, since screws are not used, focus of constraint is not possible, and that is why new fractures are a non-existent phenomenon when a flange-plate is used.

Setting the plate

The very small one or two millimeter indentions do not weaken the bone in any way. Even on senile or osteoporous bones, the result of small holes formed by the teeth is not enough to cause a focus of constraint such as one sees with the screw holes that go directly through the two cortexes and the medullary canal.

Principal advantages of this plate

The first advantage of the plate is the simplicity of its technique. Another advantage is the very low complication rate observed with this plate. Moreover, a new removal operation is not necessary for adults. In areas where patient follow-up is unsure, this is obviously an advantage. This plate is the ideal method for repairing an osteoporous bone, as well as any kind of complex fracture that often poses a problem.

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
Minimal tissue dissection for exposing the area of fracture.


When necessary, the plates may be adjusted in order to follow any kind of bone irregularity or protuberance, augmenting the space between the paired teeth by using a tooth spacer. One should not spread apart more than two pairs of teeth; instead, it's better to choose another plate.

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
Take care to conserve the periosteum during cleaning of the fractural extremities. Reduction of the fracture is done with two pliers holding the bone. They should be placed fairly far from each other, to enable the placement the Mennen plate between them.

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
Treatment of radius malunions with a specific posterior plate

The surgical treatment of radius malunions usually occurs after a long therapeutic history regularly punctuated with long periods of immobilization and professional unavailability. One must therefore propose a surgical treatment allowing to rapidly short cut the working symptomatology in order to quickly rehabilitate the functions of the hand and the wrist.

The radial fracture (lightly corticalized metaphyseal bone) is characterized by a dorso-radial compression effect, which a maladjusted osteosynthesis can perpetuate. This maladjusted osteosynthesis can even perpetuate the initial deformation and lead to a malunion. The treatment consists of an osteotomy of a postero-external opening, a cortico-spongy graft and a sufficiently stable osteosynthesis. This graft provides a basic setting, of course, but a simple synthesis with pins is not satisfactory because of a rapid osteolysis leading to creeping substitution. An anterior plate is not a good solution either because it eliminates the anterior traction forces. The posterior support and the graft are therefore in charge of absorbing the stress and strain, which is not mechanically stable. Thus, only a posterior plate permits elimination of dorsal compression, which also limits the risks of secondary movements, as much by the epiphysary grasp as by the posterior "console/bracket" effect."

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We feel that such a plate must meet very precise criteria: specific pre-molded form of minimal obstruction based upon the radius' dorsal grooves, which allows a multiple epiphyseal setting with screws or distal staples in order to minimize conflict with the extensor device, resulting in an effective stability and rapid action. This plate must be set in conformity to very strict rules with a reconstitution of the annular dorsal ligament and the reflection groove of the thumb's long extensor.

Ten patients suffered from malunions with a recurrence of the initial deformation. For five patients, the initial lesion was intra-articular and for four of them, a complementary osteotomy was needed. An anterior plate was added for a young patient because he suffered from a complex bimarginal fracture. In all cases, the cubital area was disorganized, six of which in a specific fashion (dislocation and/or fracture, TFCC lesion) requiring a follow-up gesture (osteotomy of the sigmoid cavity in de-punch fractures, arthroscopic debridement, Suvé-Kapandi operation). Good osteosynthetic stability allowed us to use eight times out of ten a lyophilized allograft rather than a cortico-spongy self-graft. The goal was actually to limit the iatrogeny of the sample area. Every time that it was possible, the dorsal annular ligament was treated between the third and the fourth area of the extensors with under-periosteal dissection in order to maintain the reflection groove of the thumb's long extensor and to protect the deep face of the screw heads' extensors. In one case, a lyophilized fascia lata patch allowed the reconstruction of a dorsal annular ligament that was damaged during previous operations. Early rehabilitation is the rule.

After at least one year, no specific complication was noted and importantly, no extensor ruptures. Only a relative weakness of the extensors could be quantified and the plates were systematically taken away. The articular criteria of reduction remained the same, outside of one case of major algodystrophy.

Concerning the surgical treatment of malunions, the radius' posterio or plate represents the most logical and stable kind of osteosynthesis. A specific pre-molded plate enables the realization of this sort of surgical intervention without any trouble, provided that the setting technique is strictly respected.

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
Fig. n° 1

Loco-regional anaesthesia

Installation on orthopaedic table under pneumatic garrot.

A dorsal arc-shaped approach with external convexity is practiced at the forearm's inferior quarter. The dorsal sensitive ramuses are carefully put aside.

One exposes the dorsal annular ligament. The latter will be opened longitudinally, just within the third sinewy area, that is to say the fourth area where one finds the common extensors.

The periosteum is incised and a centrifugal under-periosteal dissection is performed. This allows the preservation of a reflection groove of the thumb's long tendon extensor, which will be positioned anatomically during closing. During the sub-periosteal exposure, the excision of the posterior inter-osseous nerve is performed (Fig. n° 1)

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
Photo n° 2

- The bone lesions are then exposed. A preoperative tracing will allow the determination of the place and stature of the posterior opening osteotomy (Photo n° 2).

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Photo n° 3

- An autogeneous graft or a pre-cut lyophilized allograft of ideal size is then inserted with the graft-punch (Photo n° 3).

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Photo n° 4

- A per-operative X-ray is performed in order to check the correct reduction of the radial malunion and the equalization of the radio-ulnar index finger (Photo n° 4).

Photo n° 5

The pre-molded posterior plate is then positioned, starting with the upper metaphyseal oval opening, which demands a certain degree of rotation and axial sliding in order to make the implant's positioning the most appropriate possible.

Ideally, it rests on the whole radial metaphysis and should be level distally with the radius' posterior border.

The distal screws will be positioned taking into account their articular obliqueness, in order to avoid the radiocarpal passage.

Osteosynthesis can be finished with an oblique staple if the screws' hold is weak. A per-operative X-ray allows the verification of the setting's correctness.

The two dorsal and periosteal ligaments are folded back, in order to cover the osteosynthesis plate (Photo n° 5).

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
Fig. n° 6

The closing of the dorsal annular ligament is imperative, but should occur without excessive tension. Several artificial means can help the per-primam closing of the dorsal annular ligament :

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● Either a plastics, of the bipedicular flap type, by an axial incision of the internal part of the ligament.

● Or, a horizontal flap of the dorsal annular ligament can eventually be slid under the extensors, in order to protect them in relation to potential conflict with the plate or the screws (Fig. n° 6).

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● Finally, in annular ligaments' lesions, a lyophilized fascia lata patch can allow reconstruction of this physiological block (Photo n° 7).


Closing occurs with aspiratory drainage, usually with rapidly reabsorbable thread. A palmar splint is designed and will be kept for about eight days. After this immobilization, no added immobilization will be necessary, provided that the immediate mechanical holding is good. If the epiphyseal holding is more fragile, one will put in a resin handle for 6 weeks, starting rehabilitation immediately.

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

Stainless steel Pr. Goutallier plate for tibial osteotomy

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● *In frontal reaxes of genu varum arthrosics by tibial osteotomy, the internal opening osteotomy has the advantage of avoiding any vascular-nervous complication (paresis of the peroneal nerve, chamber syndrome(?), encountered in external closing osteotomies, and in rare arthrites noted in vault osteotomies of external setting.*

● *The substitution of the iliac tricortical grafts, which maintained the opening of the osteotomy, for a cement wedge (created at the beginning of the intervention) has permitted the elimination of complications due to the hold of iliac grafts (hæmatoma, suppuration, pains).*

● *The stabilization of the setting with an antero-internal support plate, based upon the shape of the tibia after the opening osteotomy, avoids all correction loss (the epiphysis above the osteotomy and the diaphysis below are each set with 3 cortical screws of 5 mm diameter).*

● *The bone regrowth of the osteotomy occurs from the outside to the inside without any local bone contribution. Progressive support is possible starting on the 45th day; complete support is authorized between the 70th and the 75th operative day.*

● *This technique seems to be very accurate since 75 % of the osteotomies solidify with the desired valgus (between 3° and 6° of valgus on the mechanic axes), and the greatest part of the errors of correction is only 1° from the ideal valgisation area. Finally, the plate is almost always tolerable and its ablation, in case of pain, is unusual.*

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Acetabular osteosynthesis is achieved with the use of plates and screws (3.5 and 4.5 mm \varnothing).

Fractures of the acetabulum are of very different types, because the size and the shape of the iliac bones vary significantly from patient to patient.

The various curved plates with equidistant holes that we offer are specifically designed for :


- Cross linking a large posterior fragment, the plate remaining parallel to the superior lip of the acetabulum.
- Application along the superior aspect of the pelvic brim, to set fractures of the two columns and those of the anterior column. Plates with the most pronounced radius are usually adequate for the majority of cases.

Principal advantages of these plates :


- The shape of the holes (equidistant).
- This shape permits the use of two different sizes of screw which can be placed obliquely in any direction. For maximum clearance of screw placement, we recommend 3.5 mm screws. This allows the setting of difficult-to-reach fragments.
- The plates can be easily contoured to the appropriate shape, using a bend-molder or two strong forceps.

Important practical note :

The plates should be properly applied to the columns after the fracture line has been reduced and secured by use of a forceps or intermittently placed screws. If the plate is too curved, or insufficiently curved, there is a risk of displacing the fracture and disturbing the reduction during the setting of the plate. According to each surgical situation, the plate should be bent in all directions, even according to its long axis, in order to give it a heliocoid shape adapted to the shape of the iliac bone.

Vue postérieure
Posterior view

Vue latérale
Lateral view

Vue antérieure
Anterior view