

# 9

## Pédiatrie • *Pediatrics*

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## *Closed site telescopic femoral nailing in child with imperfecta osteogenesis*

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*Telescopic femoral nailing according to Baily and Dubow, used since 25 years in our experience is an effective and reliable technic in particular to compensate fragility and femoral deformities due to imperfecta osteogenesis.*

*In order to avoid arthrotomy of the knee and if possible with all closed site nailing, we sought :*

- *A simple technique*
- *Minimal ancillary equipment*
- *Without increasing the financial cost*

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***Unlike the original technique, the implant should not be used for centromedullary drilling.***

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
### **Operative technique**

*The operative technique includes firstly percutaneous centromedullary drilling with scopic control. The drilling must be done starting perpendicular to the interval in order to avoid inducing an axial deviation of the nailed member.*

*The axial corrections are made by simple osteoclasy or by percutaneous puncturing. It is also possible to use a recent fracture site to obtain alignment of the diaphysis. It is important to properly tilt the neck of the femur.*

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

*The drill is taken out by a counter-incision at the buttock. It is then possible to place the female piece that will be progressively lifted without difficulty, given its lower diameter within the femur, using the guide to push it upward at the same time as the drill is pulled upward. Just before the entry of the female nail into the knee, the male part is put in place.*

*Here is a view with the placement of the female part and the drill which comes out by a buttock counter-incision. The female segment will be removed by the buttock counter-incision; the T piece will be screwed and crimped.*

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*Then, with a small punch that is adapted in the small notches that are made in the T pieces, the top of the nail is struck on the upper edge of the neck.*

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## ● First example

Here is an example of a correctly done nailing with a low osteotomy of the femur which made it possible, with the use of the fracture site, to properly realign this diaphysis, and the neck is in valgus. This is the position it should be in to avoid secondary mechanical complications.

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
## ● Second example

Here is another example of a nailing done with a closed site, this time simply using a recent fracture zone.

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## ● Esthetic benefit

Lastly, it is possible to do these nailings with only two small incisions of 1 cm, which is of course a major esthetic benefit with post-operation conditions that are much simpler.

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After the operation, immobilization is provided by either a very light pelvic-crural-leg cast, or a simple splint or even a simple dressing. It is important above all to avoid inducing any rotation anomaly, particularly a femoral retroversion.

The closed site telescopic nailing makes the post-operative situation simple and probably modifies the indications for the nailings which must be done earlier, before the onset of a significant femoral deformation.

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The ancillary equipment is extremely simple, with a long drill, longer than the femur to be nailed, and of a diameter slightly greater with respect to the implant that will be installed. The lower extremity of the drill includes threading that allows for the screwing on of the female part and can thus be introduced from bottom to top in the medullary canal. The second instrument is simply a long nail guide that has the diameter of the implant, allowing for pushing of the female segment within the femur.

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## *Telescopic D.M.M. nail for epiphysary attachment without opening the joints*

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
*The technique of placing the telescopic nail is simple and avoids opening the joints.*

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
### Material

*Telescopic nails are centromedullary nails locked into the epiphyses. They are composed of a female part into which slides a male part that terminates in an "L" (fig. 1). They are manufactured to size. The diameter and the length of the nail is calculated using telemetric X-rays.*

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Broche coudée / *Angled pin*


## Locking the female part

The nailing of the female part into the epiphysis is insured by a pin threaded in the end. Passed through the eyelet of the female part of the nail, the pin is attached with a screw tightened onto the threaded part of the pin. (fig. 2)

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## Locking of the male part

Many solutions are possibles :

- using a pin like the one mentioned before,
- using a standard pin,
- using a non absorbable surgical thread. The male part is then fastened to the soft tissue near the area of nail insertion.

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Mèche courte / *Short drill bits*

- One short and one long drill bit whose diameters correspond to the diameter of the nail used for drilling the bone

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Mèche longue / *Long drill bits*

- One drill bit whose diameter corresponds to that of the locking pin

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Viseur / *Pin guide*

- One pin guide used to evaluate the area of penetration of the locking pin. This guide must be adjusted pre-operatively by the surgeon. The point of entry of the pin is precisely determined by a preoperative X-ray test, once the female part is inserted

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



• The penetration takes place at the pre-spinal extra-articular surface after having moved aside the adipic ligaments underlying the patellar ligaments.

• Take the long drill bit. Drill the hole for the nail, preferably by hand or electric drill. Take into account at this time, in the open air, the necessary osteotomies (Fig. 3-4-5). The whole path of the nail should be completed with the drill bit guide; the nail is driven through the epiphysis inferior to the tibia, tested by X-ray from the front and the side in order to allow exact positioning beneath the cartilaginous plate.

• Pull out the drill bit and put the female nail in place, after having previously adjusted the guide. The locking nail is inserted into the hole of the rounded end of the female part. Drive the nail to the inferior end of the tibia (the opening should be lateral). For this step, use the impactor which tightens itself onto the femal part (Fig. 6).

• Nail the female part in the inferior epiphysis of the tibia (Fig. 7).

— The insertion point of the pin may be discerned with the aid of the drill guide: fasten the female part with the calibrated guide. Place the drill with the correct diameter (the diameter of the pin) in the guide, the bit should then be face-to-face with the bolted opening of the nail. Drill while verifying radiologically.

— Drive the pin through the inferior opening of the nail at the determined insertion point (just above the internal malleolus) from the interior toward the exterior. The pin is long enough to be manipulated by hand. Put it in position by hand, then drive it in with the chasse broche.

— Once the pin has breached the exterior of the skin, make a small incision in order to place the screw.

— Tighten the screw until the hook of the locking pin impacts completely in the internal malleolus and is buried under the periosteum (Fig. 7-8)

— Cut the surplus of the pin.

• Insert the male part of the nail into the female part (Fig. 9-10) and fasten it:



— either according to the form of the neighboring ligamentary pre-spinal surface by a non absorbable thread, threaded through the eyelet at the "L" end of the nail (metallic or nylon thread) (Fig. 11)

— or with the aid of a locking pin placed obliquely: In this case, upper part of the male part is pierced.

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- *Insertion area of the female nail : upper part of the great trochanter, moreso in its anterior section than posterior in order to avoid any vascular lesion of the femoral head.*
- *Perform the point of entry with a square tip.*
- *Insert the drill up to the center of the cartilage of the inferior growth of the femur. Perform as many osteotomies as necessary. Verify the passage of the drill by X-ray.*
- *Insert the female part of the nail after having previously adjusted the drill guide. For this, use the impactor which tightens itself onto the female part.*
- *Nail the female part into the inferior epiphysis of the femur.*
- *Insert the male part.*
- *Fasten the male part :
 
  - either by using non absorbable thread, to the soft tissue near the upper end of the femur (ligaments and aponeuroses neighboring the great trochanter). (Fig. 13).
  - or by using a locking pin: the pin is placed obliquely and fastened as described previously, using a screw. (Fig. 13)*

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### Extractions

- *Remove the locking screws,*
- *Pull out the locking pin(s) by pulling from the side of the hook, using the pin-punch,*
- *Remove the male part,*
- *Tighten the extractor onto the female part, this can be done easily even if the the bone has grown through the hole of the nail. Pull out the female part,*
- *When a child has grown a lot and the locking pin is included in the epiphysis, make a path with a drill in the direction of the hook. Go beyond the point of the hook and pull out the locking pin with pliers, turning it in the opposite direction of the tightening.*

### Avantages

The technique of placing the DMM nail is simple. It is unnecessary to open the joints at the level of the inferior end of the tibia and femur. The DMM nail allows passage to the knee solely in the prespinal area, therefore extra-articular.

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### *Distal derotation tibial plate for large fragments*

4 trous	6 trous
36.250.15	36.252.15
36.250.20	36.252.20
36.250.25	36.252.25

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*Plate to put the distal fragment of the left tibia in external rotation*

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*Plate to put the distal fragment of the right tibia in external rotation*

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### *Distal derotation tibial plate for small fragments*

4 trous	6 trous
36.254.15	36.256.15
36.254.20	36.256.20
36.254.25	36.256.25

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
## ● Plaque décalée *Pediatric screw-plate*

Nb de trous <i>Holes</i>	Long. en mm <i>Length mm</i>	Référence
2	56	34.201.02
3	56	34.201.13
3	72	34.201.03

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### *Pediatric blade-plate with pin guide*

Angle	100°	110°	120°
Longueur lame <i>Plate length</i>	Référence	Référence	Référence
35	37.210.35	37.211.35	37.202.35
40	37.210.40	37.211.40	37.202.40
45	37.210.45	37.211.45	37.202.45
50	37.210.50	37.211.50	37.202.50
55	37.210.55	37.211.55	37.202.55

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

## Varisation osteotomy of the proximal femur fixed with fork blade plate

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- Solid and stable implant
- **Simple**: the use of a marker pin and an image intensifier is not obligatory
- Reproducible technique: the equipment is simple and easy to set up
- **Tolerant**
- Limits errors of orientation of the equipment
- Does not fragilize the neck, particularly in the event of a false route
- Limitate the risk of growth plate lesion

The fork blade plate acts as a short blade plate, with its stability increased by a divergent cervical proximal screw in triangulation

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

The bending machine is composed of 3 parts : A part "body" a part "handle" 2 removable pawns, which come to be screwed on the part "body".

The pawns can be assembled on the right part or the left part of the "body" of the bending machine according to the size of the fork blade plate to be curved :

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


*Canulated screw 7.3 mm diameter  
proximally threaded*

- *It allows residual growth*
- *Easier to remove percutaneously*
- *Good stabilization*



Référence	Long. en mm <i>Length</i>
36.115.60	60
36.115.65	65
36.115.70	70
36.115.75	75
36.115.80	80
36.115.85	85
36.115.90	90
36.115.95	95
36.115.100	100
36.115.105	105
36.115.110	110
36.115.115	115
36.115.120	120

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## ● Simple and fast percutaneous implantation

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## ● The goal of treatment

- Epiphysis stabilization and limitation of avascular necrosis, chondrolysis, early arthrosis risk.
- Avoid hip dysmorphism due to epiphysiodesis of femoral neck.
- Allow planned controlateral hip fixation preventive and safe.


- In case of important residual growth, the screw may have to be changed

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- Journal of children's orthopaedics vol. 5 number 2, 83-88. Continued growth of the hip after fixation of slipped capital femoral epiphysis using a single cannulated screw with a proximal threading. Frédéric Sailhan, Aurélien Courvoisier, Océane Brunet, Franck Chotel, Jérôme Berard.

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Référence	Long. en mm <i>Length</i>
36.150.75	75
36.150.80	80
36.150.85	85
36.150.90	90
36.150.95	95
36.150.100	100
36.150.105	105
36.150.110	110
36.150.115	115
36.150.120	120

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