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INSTRUCTIONS FOR USE: OSTEOSYNTHESIS DEVICES FOR ADULTS AND CHILDREN

1. DESCRIPTION AND ASSEMBLY

All of our standard products have been CE  0120 approved since 1998: With respect to customised devices, our company complies with the requirements of Annex VIII of Directives 93/42/EEC and 2007/47/EC.

- **CORTICAL BONE (HA), CANCELLOUS BONE (HB) OR MALLEOLAR SCREWS, FOR ACETABULA OR ACETABULAR BASES:** These various types of screws can be self-tapping or cannulated (to enable use of a guide pin). They can be used with spacers.
 - **SLIPPED CAPITAL FEMORAL EPIPHYSIS SCREWS:** All screws are cannulated (to enable use of a guide pin). The H3P screw includes a dynamic screw, a nut and a compression screw.
 - **ALL PLATES, NAIL PLATES, SCREW PLATES, BLADE PLATES, SPADE PLATE** are fixed with standard bone screws. Screw-plates include 3 elements: a plate, a cannulated screw and a nut or a compression screw.
 - **ACETABULAR BASES:** fixed using screws and orthopaedic cement.
 - **HOOK FOR LARGE TROCHANTER:** attached with screws.
 - **NAILS AND PINS:** used in intramedullary nailing. Some types, such as Küntscher nails, are used singly to fix long bones whilst others, such as elastic nails, are used in fascicular nailing. Pins may be used as guides or references in certain procedures.
 - **PIN TIP:** A stop piece that slides on to a pin. The self-breaking tip holder is screwed on to the tip until breaking in order to attach the tip to the pin. The pin is then cut 2 or 3 mm from the stop.
 - **RETAINING WIRE:** malleable wire.
 - **CLIPS:** Various models are available: consult our catalogue.
 - **TELESCOPIC NAILS:** used in centromedullary nailing of the tibia and femur.
 - **BUTTONS:** disks drilled with 2 holes and fixed using surgical wire.
 - **LUQUE RODS:** Used to apply a correction force over the entire length of the vertebral column. The rods can be fixed with retaining wire.
 - **THORACIC PLATES:** fixed using retaining wire or slowly resorptive wire.
 - **BORELLY SLIDING CLIP SPLINTS:** the device is custom-configured during surgery through assembly of 3 types of stainless steel components: clips, straight splints and angled or straight connectors.
 - **PHALANGEAL INTERPOSITION CUPS:** implants positioned between the metatarsal head and the resected phalanx. These thin, concave disks are fixed to the articular capsule using surgical wire or pins.
- The blue coloration is reserved for titanium alloy implants.

2. MAJOR INDICATIONS

- **CLIPS, NAILS, PINS, PLATES, DEROTATION PLATES, NAIL PLATES, SCREW PLATES, BLADE PLATES, SPADE PLATES, SCREWS, SPACERS, HOOKS FOR LARGE TROCHANTER:** Fractures. Osteotomies. Arthrodesis. Osteosynthesis. Repeat procedures. Reconstruction after bone tumours. Congenital dislocation of the hip. Epiphysiodesis. Ligamentoplasty.
- **SLIPPED CAPITAL FEMORAL EPIPHYSIS SCREW:** slipped capital femoral epiphysis fixation for children and adolescent.
- **8-SHAPED PLATE:** Hemiepiphysiodeses for children and adolescent.
- **PIN TIP:** Epiphyseal and epiphyso-metaphyseal fractures.
- **TELESCOPIC NAILS:** Used in children suffering from imperfect osteogenesis, congenital curvature or congenital pseudarthrosis.
- **BUTTONS:** Used for transosseous sutures when the bone is soft.
- **RETAINING WIRE:** for retention of bone fragments or hip arthroplasties and fixation of medical devices.
- **LUQUE RODS:** Posterior arthrodesis of the vertebral column.
- **PHALANGEAL INTERPOSITION CUPS:** arthrosis, rheumatoid or post-traumatic arthritis, Hallux valgus, Hallux Regidus.
- **BORELLY SLIDING CLIP SPLINTS:** Thoracic traumatology: open or blunt trauma. Repair of simple or complex thoracic deformities. Parietectomy for tumours. Sternal applications in primary cardiac surgery, in coaptation or distraction. Sternalitis and septic pseudarthrosis. Transcostal thoracotomy.
- **THORACIC PLATES:** Thoracic deformities, serious trauma and reconstruction.
- **ACETABULAR BASES:** Acetabular reconstruction.

- The nature and component materials of the implant, indicated on the label, specify the recommended use for the implant.

3. CONTRA-INDICATIONS

- Severe or chronic infections, local or systemic.
- Sensitivity to materials.
- Bone destruction or demineralisation potentially affecting fixation of the implant.
- Severe muscular, neurological or vascular deficiencies affecting the limb in question.
- Any infection that could compromise the function of the implant.
- Overweight.

4. ADVERSE EFFECTS AND POSSIBLE COMPLICATIONS

The patient must be warned of the limitations and risks inherent to the placing of the implant. Some complications can require revision surgery.

- Early or late-onset infections, haematoma, allergies, cutaneous necrosis, necrosis of bone or bone resorption, neurological complications (clinical or electrical monitoring is advised: wake-up test), pain, bone resorption or fractures, mechanical complications (unscrewing, implant failure or migration, wear, loosening, etc.), implant fracture.
- Cardiovascular disorders, thrombosis, consolidation delays, pseudarthrosis, chondrolysis, damage to the growth plate (for paediatric implants), algodystrophy.

5. GENERAL INSTRUCTIONS FOR USE

- Pre-operative planning:

STANDARD IMPLANTS: The choice of the type and characteristics (angle, dimensions) of the implant to be used depends on patient pathology and morphology.

CUSTOMISED IMPLANTS: Pre-operative planning has provided the essential information required for manufacture of a customised implant. The device was designed for a patient with help from the surgeon. Instructions concerning pre-operative planning and overlays are provided with the implant.

- Surgical technique: see sales documentation. Implants should be placed using the ancillary material provided for this purpose. Please contact us regarding customised implants.

 : Attention

- For more information on the characteristics of the drills, taps and cutters to be used with screws, please see standard NF S 90-432 or ISO 9714-1.
- Each customised implant is designed for a specific patient and should never be used on another patient. Ensure that the implant corresponds to the indications appearing on the packaging and overlay.
- Don't use a titanium screw with a stainless steel implant (nor vice versa), as this would accelerate the process of corrosion.
- **BORRELLY SLIDING CLIP SPLINTS:** the hinges should be crimped very carefully to minimise any risk of decrimping and implant migration. Crimping should be performed before modifying the faces or sides of the curved pieces in any way. If the hinge at the clip or connector is not sliding optimally on the splint then it can be enlarged slightly using decrimping forceps. It is important to mark at least two crimping points for each two-component assembly. The decrimping forceps for the hinges should be properly sharpened. Return them to us for repair if necessary.
- Perform a visual inspection of the implant to detect any possible damage (scratches, pitting, etc.). Do not allow the implant to come into contact with instruments that could mark its surface. Do not use an implant if it has been subject to any impact or has visible damage.
- Our company cannot be held liable for any corrective modifications made to implants.
- Improper positioning of an implant could result in its failure. Implants should therefore be placed only by a surgeon with the necessary training.
- Do not use implants as test devices.
- Dispose of implants and waste coming into contact with patients in accordance with current regulations.

Some implants are made of multiple components. Do not combine components from different manufacturers or origins. In this event, our company waives all liability.

: Do not use if the packaging is damaged

: Do not reuse

6. POST-OPERATORY PRECAUTIONS

Advise the patient of the precautions to be followed as part of post-operative development.

- Carry out radiological surveillance (frequency and protocol determined by surgeon).
 - Proceed with timely and effective treatment of any infection (nearby cutaneous, urinary, bronchopulmonary, dental, etc.), however benign, due to risk of haematogenous contamination.
 - Be attentive to any signs of pain.
 - Consult us if revision surgery appears necessary.
- The surgeon is to decide on authorising free movement of the patient and to define any limits on patient activity after implantation. However, excessive activity involving the affected limb is strongly discouraged.

7. ABLATION OF DEVICE

The surgeon makes the final decision concerning ablation of an implant.

8. INTERFERENCE WITH OTHER MEDICAL TESTS OR TREATMENTS

- The presence of stainless steel in the implant generates artefacts in the immediate vicinity of the device. This must be taken into account in relation to medical imaging tests.
- Intramuscular injections on the same side as the affected limb are strongly discouraged.

9. HANDLING-STORAGE

Follow the required precautions for asepsis instructions when removing from packaging.

Some products have sharp components that could injure the handler.

The packaging of **products delivered sterile** is not to be opened until time of use. Store products in their packaging.

: fragile, handle with care.

10. STERILITY

- **Products delivered sterile:**

: Sterilized by irradiation

Re-sterilization of products delivered sterile by our company may cause risks of infection and/or cross-contamination.

- **Products delivered non-sterile:**

: Non-sterile

Non-sterile implants are to be sterilized within the health care facility. They are delivered in a protective bag which cannot be used as sterile packaging: they must be cleaned, decontaminated and sterilized according to the regulations in effect. We recommend a sterilization by autoclave at 134°C during minimum 18 minutes.

Repeated moist heat sterilization of implants have minimal effects on implants. However, implants must systematically be inspected. Any implant with corrosion, scratches, notches, residue, discoloration or debris should be discarded.

11. MEANING OF THE SYMBOLS

REF	Catalogue number		Manufacturer		Consult precautions for use
LOT	Batch number		Date of manufacture		Limit date for use

For additional information, please contact us or your dealer or consult the videos and documentation available at our website www.medicalex.info

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