1. DESCRIPTION

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1. DESCRIPTION

Cranial halo : It is composed of an halo hoop fixed to the skull with 4 or 5 pins. In particular case of a very young child, the surgeon could use 8 pins. It can be combined with traction cradles or a thoracic contention device.

Cardner-Wells tractor : It is made of a traction arc pinned to the skull with two pins.

External braces : They have a certain number of standards parts: pins, some joints, some ball-and-socket hafts, some union and sliding rods. The mobility of the rods inside the ball-and-socket joint enables adaptation of the brace to most anatomical situations.

For the Pins and the wires :  : Do not reuse.

Other elements are reusable. The constituent material is indicated on the label.

CE labelling since 1998 :


*Non implantable components:  * Pins and wires : 

For custom-made products, our company conforms to the requirements of the additional clause VIII of the Directive 93/42/CEE.

2. ASSEMBLY

Do not assemble elements of devices from diverse origins. Use only components produced by our company. If not, we can not be held responsible. Titanium pieces are not compatible with the stainless steel pieces and conversely.

Cardner-Wells tractor's assembling drawing :

3. MAIN INDICATIONS

Halo and Tractor : Spinal fractures and luxations, fusions, spinal deformations, scolioses, tumors, rheumatoid arthritis.

Braces : Fractures, loss of bone substance, arthrodeses, osteotomies, persistent calluses, pseudoarthroses, elongations, osseous malformations, resumption of previous interventions.

4. COUNTER-INDICATIONS

Genito-urinary pathology, occipito-atlantoidienne luxation, osseous mineral deficiency which could affect the device fixing, disease of the skin, sensitivity to the material.

Severe neurological, vascular or muscular deficiencies, which affect the concerned extremity.

5. SIDE-EFFECTS AND POSSIBLE COMPLICATIONS

The patient should be informed of the limits and danger inherent in the placing of the devices mentioned above. Some complications may lead to a change in the pins, wires, or to ablation of the device.

- Pains, rigidity,
- Neurological problems,
- Cutaneous problems (scabs, necroses, post-operative scars), allergies, internal or external infections,
- Bone lysis, pseudoarthroses, slowed consolidation, persistent calluses,
- Mechanical complications : material ruptures, dismantling,
- Thrombosis, cardio-vascular disorders,
- Deformation of the skull, penetration of the pins at the internal table, vertebral subluxations, intracranial hypertension, desarthrodeses.

6. INSTRUCTIONS FOR GENERAL USE

- **Operative planning**: The choice of the device to implant must be made according to the morphology and pathology of the patient.

- **Operative technique**: see sales documentation.

The device must be set with the ancillary material provided for this.

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Cardner-Wells Tractor : The hair and the pins have to be cleaned with an anti-septic aerosol. The penetration point of the pins is marked on the vertical axis of the external auditory canal, at the temporal crest. At each penetration point, local anesthesia is injected. After a few moments, the tractor is applied with or without skin incision. Then, the system is retightened until the pressure indicator projects to 1mm. Then, the tractor is moved back and forth to stabilize the pins in the skull.

: Attention

- Check that the product corresponds to the indications on the packaging.
- Follow required asepsis directions when removing from packaging.
- Proceed with a visual exam of the pins and wires in order to detect potential deteriorations. In case of visible defects (i.e.: bent pin points or wires), implants must not be used.
- In case of corrective actions on devices, our company can not be held responsible.
- Do not apply the pins on fragile parts of the skull.
- The device must be set by a surgeon who has acquired the necessary training.
- Do not ever re-use pins or wires which have been already used.

7. POSTOPERATIVE PRECAUTIONS

- Perform a radiological surveillance whose periodicity and protocol will be defined by the surgeon.
- Perform a neurological surveillance. (eyes mobility, swallowing, lingual movements)
- Perform a local surveillance at the penetration points of the pins and/or of the wires.
- Check regularly that the pins are well-tightened,
- Perform an orthopaedic surveillance, especially for rigid fixations,
- In cases of halo-casts, when the patient is lying down, the halo must not lean on the bed.

It's up to the surgeon to define the patient's activity.

8. ABLATION OF THE MATERIAL

It's up to the surgeon to make the final decision in regards to the ablation of the device.

9. INTERFERENCES WITH OTHER EXAMS OR MEDICAL TREATMENTS

The presence of stainless steel generates artefacts around the material. Take it in consideration during medical imaging exams.

9. HANDLING-STORAGE

Some products have sharp components that could injure the handler.

: Fragile, handle with care.

10. STERILITY

All the material is sold not sterile, with the exception of the pins which can be sold sterile or not sterile.

- Non-sterile devices are to be sterilised within the health care facility. They are delivered in a protective bag which cannot be used as sterile packaging: they must be cleaned, decontaminated and sterilised according to the regulations in effect. We recommend a sterilisation by autoclave at 134°C during minimum 20 minutes.
- For the pins sterilised :

: Do not re-sterilize.  : Sterilized by irradiation.

: Do not use if the packaging is damaged.

11. MEANING OF THE SYMBOLS

 REF	Catalog reference	 Manufacturer	 Consult precautions for use
 LOT	Batch's number	 Date of manufacture	 Limit date for use

For additional information, please contact us or your dealer or consult the videos and documentation available at our website www.medicalex.info