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INSTRUCTIONS FOR ALL ANCILLARY, INSTRUMENTS AND REUSABLE MATERIALS

DECONTAMINATION-CLEANING-STERILIZATION

- Non-sterile devices are to be sterilized within the healthcare facility. They are delivered in a protective bag which cannot be used as sterile packaging. They must be cleaned, decontaminated and sterilized according to the regulations in effect. We recommend a sterilization by autoclave at 134°C during minimum 20 minutes.
- Respect the asepsis instructions found during unpacking.
- In case of lending, products must be sent back decontaminated and cleaned.
- For efficient stages of cleaning/decontamination, do not forget to:
 - ✓ Dismantle the products which can be. (In case of doubt, contact us)
 - ✓ Clean with a swab inside of the drilled areas.
 - ✓ Brush the parts of the products which need it (screw threading, cavity, etc.)
- Dispose of waste and implants that have come into contact with patients according to current regulations.

MATERIAL USED WITH A MOTOR

- Use a motor dedicated to surgery.
- Check the tightening of our device in the mandril before turning on the motor.
- Refer to the instructions of the motor in order to learn how it functions before using it in association with our device.

DEVICES WITH A MEASURING FONCTION

Limits of accuracy of the Weil measuring forceps : $\pm 0.25\text{mm}$

NOTE REGARDING PEDIATRIC BENDING MACHINES

Instructions for assembly and disassembly are available in the catalog and on the website www.medicalex.info
 Graduations on the bending machines are given as a rough guide and do not provide accurate measurements.
 The bending machine for the fork plate (Ref 35.600.04) must only be used with fork plates from our company.
 The bending machine for the H3P screw-plate (Ref 34.330.60) must only be used with H3P screw- plates from our company.

WARNINGS

Attention

- Surgeons and staff using the devices must have received the necessary training.
- Perform a visual inspection of the device in order to detect any possible damage (scratches, pitting...). Do not use a device if it has been subjected to any impact or has visible damage.
- Do not allow the device to come into contact with surfaces or instruments that could mark its surface.
- As soon as an anomaly is noted, please contact the company in order to take suitable measures.

HANDLING

Fragile, handle with care.

Some products are sharp and can injure the handlers.

MEANING OF THE SYMBOLS

REF	Catalog reference		Consult precautions for use	LOT	Batch's number		Manufacturer		Non sterile
------------	----------------------	---	-----------------------------------	------------	-------------------	--	--------------	---	----------------

For additional information, especially regarding surgical techniques, please consult the videos and documentation available at our website www.medicalex.info